

# ANL SERIES

## High Energy and High Repetition Rate DPSS Nanosecond Lasers

Highly Customizable  
to Meet Customer  
Needs



ANL series electro-optically Q-switched nanosecond Nd:YAG lasers deliver high energy pulses at high repetition rates.

A diode-pumped Q-switched nanosecond laser, based on industry-tested technology is used as a master oscillator of the system. It produces high-intensity, high-brightness pulses and is well suited for further amplification in linear amplifiers for high-energy flat-top output pulses. Employing electro-optical cavity dumping, the master oscillator can produce pulses which are as short as several ns with uniform beam profile and low divergence.

Power amplifiers are a chain of low-maintenance diode-pumped single and double pass amplifiers

where pulses are amplified up to the required energy. During amplification, spatial beam shaping is employed in order to get a flat top shape at the output. Optional second and third harmonic generators are based on angle-tuned nonlinear crystals placed in heaters.

For convenience, PC software for Windows™ (LabVIEW™ drivers are supplied as well) is used for laser operation, monitoring and internal system diagnostics.

To tailor the laser for specific applications or requirements, various customization possibilities are available such as industrial grade, portable laser housing with integrated power supplies and cooling units.


### FEATURES

- ▶ Up to **1 J** at **1064 nm** output pulse energy
- ▶ Up to **1 kHz** repetition rate
- ▶ **2 – 4 ns** or **5 ns** pulse duration
- ▶ Spatial flat top beam profile
- ▶ Low maintenance costs
- ▶ Various customizing possibilities to tailor for specific applications
- ▶ Optional second and third harmonics generators
- ▶ High efficiency diode pumping chambers
- ▶ **1×2 m** laser head footprint
- ▶ Internal system diagnostics
- ▶ Optional industrial grade, portable laser housing with integrated power supplies and cooling units

**SPECIFICATIONS <sup>1)</sup>**

Model	ANL 2001k	ANL 4001k	ANL 1k200
<b>MAIN SPECIFICATIONS</b>			
Pulse energy			
at 1064 nm	> 200 mJ	> 400 mJ	> 1000 mJ
at 532 nm <sup>2)</sup>		–	
Pulse energy stability (StdDev): <sup>3)</sup>			
at 1064 nm		1.5 %	
at 532 nm		–	
Power drift <sup>4)</sup>		± 2 %	
Pulse duration <sup>5)</sup>		2 – 4 ns	~ 5 ns
Repetition rate		1000 Hz	200 Hz
Polarization at 1064 nm		horizontal	
Optical pulse jitter <sup>6)</sup>		–	
Linewidth		–	
Beam profile		Hat-Top (at laser output), without diffraction rings	
Typical beam diameter <sup>7)</sup>		~6 mm	~10 mm
Beam divergence <sup>8)</sup>		< 1.0 mrad	< 0.5 mrad
Beam pointing stability		± 30 µrad <sup>3)</sup>	
<b>PHYSICAL CHARACTERISTICS</b>			
Laser head (W × L × H)		1000 × 2000 × 490 mm	
Power supply unit (W × L × H)		553 × 600 × 700 mm	
Umbilical length		2.5 m	
<b>OPERATING REQUIREMENTS</b>			
Facility water consumption (max 20° C)	10 l/min	14 l/min	10 l/min
Ambient temperature		22 ± 2 °C	
Relative humidity		20 – 80 % (non-condensing)	
Power requirements <sup>9)</sup>		208, 380 or 400 V AC, three phase, 50/60 Hz	
Power consumption	<10 kW	<12 kW	<6 kW

- <sup>1)</sup> Due to continuous improvement, all specifications subject to change without notice. Parameters marked typical may vary with each unit we manufacture. Unless stated otherwise, all specifications are measured at 1064 nm and for basic system without options.
- <sup>2)</sup> For NL94X-SH harmonic generator option. Harmonic outputs are not simultaneous; only single wavelength beam is present at the output at once.
- <sup>3)</sup> Standard deviation value averaged over 30 s after 20 minutes of warm-up.
- <sup>4)</sup> Deviation from average value measured over 8 hours of operation when room temperature variation is less than ±2 °C.
- <sup>5)</sup> Measured with photodiode with 100 ps rise time and oscilloscope with 600 MHz bandwidth.
- <sup>6)</sup> Standard deviation value, measured with respect to triggering pulse.
- <sup>7)</sup> Beam diameter is measured at 1064 nm at laser output at the 1/e<sup>2</sup> level and can vary with each unit we manufacture.
- <sup>8)</sup> Full angle measured at the 1/e<sup>2</sup> level at 1064 nm.
- <sup>9)</sup> Mains voltage should be specified when ordering.



Typical beam profile of ANL4001k laser




[www.phototechnica.co.jp](http://www.phototechnica.co.jp)  
**フォトテクニカ株式会社**  
 〒336-0017 埼玉県さいたま市南区南浦和 1-2-17  
 TEL:048-871-0067 FAX:048-871-0068  
 e-mail:voc@phototechnica.co.jp

# PGx01 SERIES



## High Energy Broadly Tunable OPA

### FEATURES

- ▶ Ultra-wide spectral range from **193 to 16000 nm**
- ▶ High peak power (>**50 MW**) ideal for non-linear spectroscopy applications
- ▶ Narrow linewidth <**6 cm<sup>-1</sup>** (for UV < 9 cm<sup>-1</sup>)
- ▶ Motorized hands-free tuning in 193–2300 nm or 2300–16000 nm range
- ▶ PC control via USB port (RS232 is optional) and LabVIEW™ drivers
- ▶ Remote control via keypad

Travelling Wave Optical Parametric Generators (TWOPG) are an excellent choice for researchers who need an ultra-fast tunable coherent light source from UV to mid IR.

### Design

The units can be divided into several functional modules:

- ▶ optical parametric generator (OPG);
- ▶ diffraction grating based linewidth narrowing system (LNS);
- ▶ optical parametric amplifier (OPA);
- ▶ electronic control unit.

The purpose of the OPG module is to generate parametric superfluorescence (PS). Spectral properties of the PS are determined by the properties of a nonlinear crystal and usually vary with the generated wavelength. In order to produce narrowband radiation, the output from OPG is narrowed by LNS down to 6 cm<sup>-1</sup> and then used to seed OPA.

Output wavelength tuning is achieved by changing the angle of the nonlinear crystal(s) and grating. To ensure exceptional wavelength reproducibility, computerized control unit driven precise stepper motors rotate the nonlinear crystals and

diffraction grating. Nonlinear crystal temperature stabilization ensures long-term stability of the output radiation wavelength.

In order to protect nonlinear crystals from damage, the pump pulse energy is monitored by built-in photodetectors, and the control unit produces an alert signal when pump pulse energy exceeds the preset value.

For customer convenience the laser can be operated from master device or personal computer through USB (VCP, ASCII commands), RS232 (ASCII commands) or LAN (REST API) interfaces or from remote control pad with backlit display that is easy to read even while wearing laser safety glasses.

### APPLICATIONS

- ▶ Nonlinear spectroscopy: vibrational-SFG, surface-SH, Z-scan
- ▶ Pump-probe experiments
- ▶ Laser-induced fluorescence (LIF)
- ▶ Other laser spectroscopy applications

### Available models

Model	Features
PG401	Model has a tuning range from 420 to 2300 nm and is optimized for providing highest pulse energy in the visible part of the spectrum. The wide tuning range makes PG401 units suitable for many spectroscopy application.
PG501-DFG	Model has a tuning range from 2300 to 16000 nm. The PG501-DFG1 model is the optimal choice for vibrational-SFG spectroscopy setups.

**SPECIFICATIONS <sup>1)</sup>**

Model	PG401	PG401-SH	PG401-DUV	PG501-DFG1	PG501-DFG2
<b>Tuning range</b>					
DUV	-		193–209.95 nm	-	
SH	-	210–340, 370–419 nm	-		
Signal	420 – 680 nm	-			
Idler	740 – 2300 nm	-			
DFG				2300–10000 nm	2300–16000 nm
Output pulse energy <sup>2)</sup>	> 1000 µJ at 450 nm	> 100 µJ at 300 nm	> 50 µJ at 200 nm	> 250 µJ at 3700 nm, > 40 µJ at 10000 nm	> 250 µJ at 3700 nm, > 80 µJ at 10000 nm
Linewidth	< 6 cm <sup>-1</sup>	< 9 cm <sup>-1</sup>		< 6 cm <sup>-1</sup>	
Max pulse repetition rate	50 Hz				
<b>Scanning step</b>					
Signal	0.1 nm	-			
Idler	1 nm	-			
Typical beam size <sup>3)</sup>	~4 mm	~3 mm		~9 mm	
Beam divergence <sup>4)</sup>	< 2 mrad			-	
Beam polarization	-	vertical		horizontal	
Signal	horizontal	-			
Idler	horizontal	-			
Typical pulse duration	~15 ps	~12 ps		~20 ps	

**PUMP LASER REQUIREMENTS**

<b>Pump energy</b>					
at 355 nm	-	10 mJ		-	
at 532 nm	-			10 mJ	
at 1064 nm	-	2 mJ	6 mJ	15 mJ	
Recommended pump source <sup>5)</sup>	PL2231-50-TH, PL2251A-TH		PL2231-50-TH, PL2251A-TH	PL2231A-50-SH, PL2251B-SH	
Beam divergence	< 0.5 mrad				
Beam profile	homogeneous, without hot spots, Gaussian fit >90 %				
Pulse duration <sup>6)</sup>	30 ± 5 ps				

**PHYSICAL CHARACTERISTICS**

Size (W x L x H)	456 × 633 × 244 mm	456 × 1031 × 249 ± 3 mm			
------------------	-----------------------	-------------------------	--	--	--

**OPERATING REQUIREMENTS**

Room temperature	15 – 30 °C				
Power requirements	100 – 240 V AC single phase, 47 – 63 Hz				
Power consumption	< 100 W				

<sup>1)</sup> Due to continuous improvement, all specifications are subject to change without notice. Parameters marked typical are not specifications. They are indications of typical performance and will vary with each unit we manufacture. Unless stated otherwise, all specifications are measured at 450 nm for PG401 units, 3000 nm for PG501 units and 300 nm for PG401SH units and for basic system without options.

<sup>2)</sup> See tuning curves for typical pulse energies at other wavelengths. Higher energies are available, please contact Ekspla for more details.

<sup>3)</sup> Beam diameter is measured at the 1/e<sup>2</sup> level.

<sup>4)</sup> Full angle measured at the FWHM point.

<sup>5)</sup> If a pump laser other than PL2250 or PL2230 is used, measured beam profile data should be presented when ordering.

<sup>6)</sup> Should be specified if non-EKSPLA pump laser is used.



CUSTOMIZED FOR SPECIFIC REQUIREMENTS

Please note that these products are custom solutions tailored for specific applications or specific requirements.

Interested? Tell us more about your needs and we will be happy to provide you with tailored solution.

PG401-DFG1 provides:

- ▶ The broadest hands-free tuning range – from 420 to 10000 nm
- ▶ It can be further extended up to 16000 nm with -DFG2 option. It should be noted, that for the 8000 – 16000 nm range a different nonlinear crystal is used, and exchange of the crystals needs to be done manually

PG402 features:

- ▶ Gap-free tuning range 410 – 709, 710 – 2300 nm
- ▶ Linewidth < 18 cm<sup>-1</sup>

TUNING CURVES




Fig 1. Typical PG401 model tuning curve  
Pump energy: 10 mJ at 355 nm




Fig 2. Typical PG401-DUV model tuning curve




Fig 3. Typical PG401-SH model tuning curve. Pump energy: 10 mJ at 355 nm

Note: The energy tuning curves are affected by air absorption due narrow linewidth. These pictures present pulse energies where air absorption is negligible.




Fig 4. Typical PG501-DFG1 tuning curve in 2300–10000 nm range  
Pump energy: 7 mJ at 1064 nm




Fig 5. Typical PG501-DFG2 tuning curve in 8000–16000 nm range  
Pump energy: 15 mJ at 1064 nm

RECOMMENDED UNITS ARRANGEMENT ON OPTICAL TABLE




Fig 6. Arrangement of pump laser and PGx01 unit on optical table




Fig 7. Recommended arrangement of pump laser and PGx01-DFGx unit on optical table

OUTLINE DRAWINGS




Fig 8. PG401 external dimensions

OUTPUTS PORTS

Model	L, mm	a, mm	b, mm	c, mm	Port 1	Port 2
PG401	633	380	x	x	420–680 nm, 740–2300 nm	–
PG401-SH	838	380	x	x	210–340 nm, 370–419.9 nm, 420–680 nm, 740–2300 nm	–
PG401-SH/DUV	1026	380	250	50	210–340 nm, 370–419 nm, 420–680 nm, 740–2300 nm	192–209.95 nm




Fig 9. PG501 external dimensions

For SFG optional 532 nm output port 2.

ORDERING INFORMATION

**PG401-DUV**

Model  
PG4xx → 355 nm pump

01 → travelling wave, narrowed linewidth  
02 → travelling wave, not narrowed  
11 → synchronous pumping, narrowed

Optional tuning range extension  
DUV → 193–209.95 nm  
SH → 210–340 nm & 370–420 nm

*Custom products, tailored for specific applications. Inquire for other specifications.*

DFG1 → 2300–10000 nm; >250 μJ at 3700 nm  
DFG2 → 2300–16000 nm

**PG501-DFG1**

Model  
PG5xx → 532 nm pump

01 → travelling wave, narrowed linewidth

Tuning range  
DFG1 → 2300–10000 nm; >250 μJ at 3700 nm  
DFG2 → 2300–16000 nm

**Note:** Laser must be connected to the mains electricity all the time. If there will be no mains electricity for longer than 1 hour then laser (system) needs warm up for a few hours before switching on.

**PHOTO TECHNICA** フォトテクニカ株式会社  
www.phototechnica.co.jp  
〒336-0017 埼玉県さいたま市南区南浦和 1-2-17  
TEL:048-871-0067 FAX:048-871-0068  
e-mail:voc@phototechnica.co.jp